



## EQUINE DYSTOCIA SIMULATOR

WITH EQUINE NASOGASTRIC INTUBATION AND NECK VENIPUNCTURE  
W/ INTRAMUSCULAR INJECTION

The VSI Equine Dystocia Simulator is the latest addition to our suite of veterinary simulators aimed at improving clinical outcomes. Designed with accuracy and functionality in mind, both mare and foal have been modelled to be realistic in size and conformation. The foal is equipped with an articulated skeleton which allows it to be placed in a variety of positions allowing students to gain the confidence needed to perform dystocia scenarios in a live setting. Nasogastric intubation and venipuncture simulations can also be performed as these features have been added to our latest model. The VSI equine uterus set is compatible with the VSI Equine Dystocia Simulator resulting in a versatile and practical addition to any equine clinical skills training program.



## EQUINE BODY FEATURES

- Modeled as a standard 15 hand Quarter Horse
- Steel reinforced epoxy/Fiberglass construction with water resistant components for ease of cleaning
- 1.36m at the shoulder and 2.44m from nose to tail. Widest point is 0.8m
- Inflatable foal air bed support system
- Clear vinyl uterine bag
- Soft durable birthing perineum panel
- Replica pelvis
- Padded fetal extractor, obstetric chain, head snare and landing mat
- Abdominocentesis function with replaceable pad; fluid reservoir and palpable linea

## REPLACEABLE PARTS

- Birthing Panel Assembly 1025591
- Air Bed 1025588
- Pelvic Air Bladder 1025598
- Foal 1025587



## DAY ONE SKILLS & PRODUCT BENEFITS

Students will gain confidence and develop familiarity with foal positioning and palpation.

- 34kg (75lbs) birthing foal size
- Airway for intubation tube
- Rectal palpation
- Durable, flexible and easy to clean skin
- Compatible with the equine uterus set for artificial insemination training upgrade kit and palpation of multi stage ovaries
- Realistic in size and aesthetically pleasing allowing for better immersion in training scenarios

The Equine Nasogastric Intubation and Neck Venipuncture/IM Injection is removable from the equine body and may be used as a stand alone unit. Can be purchased separately.

### Features include:

- Jugular venipuncture with palpable veins that can endure multiple punctures (20 or 22 gauge needle recommended for simulation)
- Easily replaceable latex vein tubing that can endure multiple punctures; includes 5 sets of veins
- Mechanical pressurization system allows for positive pressure for blood draw
- Long lasting, easily removable, washable, replaceable hide covering with mane
- IM site that can endure multiple punctures and accept fluid injections
- Soft silicone head with flexible ears and well rendered anatomy illustrates facial injection sites and accepts a bridle
- Rolling Stand, adjustable for height, rostral tilt, and 90 degree rotation.

**Equine Dystocia Simulator 1025586**

## UPGRADABLE PARTS

Artificial Insemination techniques can be practiced with the purchase of the Palpation Panel With Equine Theriogenology Uterus Set. Equine uterus set includes interchangeable ovaries.



IF YOU HAVE ANY QUESTIONS OR NEED SUPPORT, PLEASE SCAN HERE



**3B Scientific Group**

Headquarters · Ludwig-Erhard-Strasse 20 · 20459 Hamburg · Germany · [3bscientific.com](http://3bscientific.com)