



**STOPS
MEDICAL**
A 3B SCIENTIFIC COMPANY

3B
3B Scientific



6-IN-1 TCCC TRAINER


The 6-in-1 Tactical Combat Casualty Care (TCCC) Trainer is comprised of a Hyper-Realistic®, anatomically similar male head and upper torso. It is constructed with a skeleton that utilizes user-repairable body skin, neck skins, tracheas, and user-replaceable intraosseous training pucks.

Procedures available on the 6-in-1 TCCC Trainer include the following:


1. Nasopharyngeal Airway (NPA) Adjunct Insertion (with Oropharyngeal visual confirmation of proper measurement and placement)
2. Clearing upper airway of obstructions (i.e. finger-sweep) with Oropharyngeal airway measurement and placement
3. Surgical airways (Cricothyroidotomy, Cricothyroidcentesis, etc.)
4. Chest Needle Thoracentesis, 4 Sites, MCL, AAL
5. Sternal Intraosseous Catheter Insertion
6. Proximal Humerus Intraosseous Catheter Insertion




Nasopharyngeal Airway (NPA) Adjunct



Clearing upper airway of obstructions




Surgical airways



Chest Needle Thoracentesis



Sternal Intraosseous Catheter Insertion



Proximal Humerus Intraosseous Catheter Insertion

6-IN-1 TCCC Trainer 1025848

IF YOU HAVE ANY QUESTIONS OR NEED SUPPORT, PLEASE SCAN HERE



02/2025



3B Scientific Group

Headquarters · Ludwig-Erhard-Strasse 20 · 20459 Hamburg · Germany · 3bscientific.com